

ROTORS				ROTORS			
10 mm XC10 Length =35.00 mm				10 mm SuperSonic Length =27.50 mm			
Sample volume – .6 mL to 1.10 mL				Sample volume – .6 mL to 1.10 mL			
Low density* Spinning maximum kHz				Low density* Spinning maximum kHz			
# 44265	Zirconia XC Thick Wall Rotor	8.5	\$735	# 42113	Silicon Nitride SS Thick Wall Rotor	8.5	\$ 725
# 44266	Zirconia XC Thin Wall Rotor	4.5	775	# 42138	Zirconia SS Thick Wall Rotor	6	725
				# 42193	Silicon Nitride SS Thin Wall Rotor	8	775
				# 42173	Zirconia SS Thin Wall Rotor	4.5	775

10 mm XC/SS Part #	Description	Price
Caps for THICK Wall Rotors – fastest spinning		
96436	XC/SS Kel-F Short Thick Cap Pair	\$155
-	XC Only Kel-F Long Thick Cap Pair	-
-	XC/SS GFT Short Thick Cap Pair (glass filled torlon)	-
-	XC Only GFT Long Thick Cap Pair (glass filled torlon)	-
96438	XC/SS Aurum or Torlon Short Thick Cap Pair	175
-	XC Only Aurum or Torlon Long Thick Cap Pair	-
96449	Kel-F O-Ring Cap Pair	250
96451	Aurum or Torlon O-Ring Cap Pair	250
Caps for THIN Wall Rotors		
96471	XC/SS Kel-F Short Thin Cap Pair	\$155
-	XC Only Kel-F Long Thin Cap Pair	-
-	XC/SS GFT Short Thin Cap Pair (glass filled torlon)	-
-	XC Only GFT Long Thin Cap Pair (glass filled torlon)	-
96473	XC/SS Aurum or Torlon Short Thin Cap Pair	175
-	XC Only Aurum or Torlon Long Thin Cap Pair	-
Cap Pullers and Accessories		
96174	Puller for Turbine Caps	\$210
96190	Threaded Cap Puller for O-ring Caps	110
-	Rotor Holder – use with plungers below	-
-	Plunger - thick wall- to push out caps or sealing cells	-
-	Plunger - thin wall- to push out caps or sealing cells	-
96508	XC Sample Packing Set for Thick Wall Rotors	\$160
96507	XC Sample Packing Set for Thin Wall Rotors	160
96508	SS Sample Packing Set for Thick Wall Rotors	\$160
96507	SS Sample Packing Set for Thin Wall Rotors	160